

Approved Manufacturers

Clinell Peracetic Acid Wipes* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Abecca	Positioning Supports and Mattresses
Active Key GmbH & Co. KG	All Active Key MedicalKey™ products, including: AK-C4112, AK-C4412 AK-C7012 and AK-C7412
AGFA	DX-D products
Agile Medical	All products
BD (Bard)	Site~Rite® 8 Ultrasound System and Probe
Carrflex	Mattress Covers, Upholstery Tops and Associated Equipment
Carr Plastics	Techmaflex
Dartex Coating	Mattresses (with polyurethane face)
Fresenius Kabi - Agilia range	VP MC Volumetric Infusion Pump SP MC Syringe Infusion Pump SP TIVA TCI Syringe Infusion Pump SP PCA Syringe Infusion Pump
Fujifilm Healthcare	Ultrasound Probes and Accessories (suitable for most, please see cleaning instruction QAD-Z3510 for full details)
Fujifilm Sonosite	Ultrasound Equipment (Edge, NanoMaxx, MicroMaxx, TITAN) iViz Devices and C60v, L25v, L38v and P21v probe
Healthy Mind	Pico 4 Enterprise Pico G2 4K Enterprise Associated Tablet Lenovo M10 Associated Bose headphones
Howard Wright	Mattresses

*Previously Clinell Sporidical Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 16.

www.gamahealthcare.com/compatibility

Approved Manufacturers

Clinell Peracetic Acid Wipes* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
HP	HP EliteBook 840 G6 Healthcare Edition HP EliteOne 800 G5 Healthcare Edition All-in-One HP Healthcare Edition HC270cr Clinical Review Display HP Healthcare Edition HC241 Clinical Review Monitor HP Healthcare Edition HC241p Clinical Review Monitor HP Healthcare Edition HC271 Clinical Review Monitor HP Healthcare Edition HC271p Clinical Review Monitor
Mindray	N Series Monitor, HyLED X Series, HyLED C Series, HyPort P Series, HyPort R80 Series and HyBase V Series
Modsel	All product ranges
OES Medical	Pressure Monitor and UAM Anaesthetic Machine
Orfit Industries	All related types of AIO 3.0 Medium Density Blue Head Support Low Density Blue Head Support
Panaz	Cadet Furniture
Parity Medical	Mobile Computer Carts

*Previously Clinell Sporocidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 16.

www.gamahealthcare.com/compatibility

Approved Manufacturers

Clinell Peracetic Acid Wipes* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Pegasus Medical	Poly Line (Cart Types: Treatment, Emergency, Drug Distribution, Anaesthesia & A.I.O)
	Ecoline Line (Treatment, Emergency, Drug Distribution, Anaesthesia & A.I.O)
	Apollo Line (Catheter, Medication, General Storage, MRI, Sutures)
	MaxSalus Line (All)
	TopSalus Line (All)
	Hygieia Line (SPD/CSSD)
	Steel Rack Line (All)
	Roll Stand Lines (Medical Devices, Telemedia, IT, Blood Samples)
Philips	ISO 4060 Storage Modules (Antibacterial Module Line, ABS Module Line, Polycarbonate Module Line, High Resistance Module Line, Stainless Steel Module Line)
	Ultrasound Transducers (please email info@gamahealthcare.com for full list of compatible transducers)
Recticel	Mattresses
Sidhil	IQ Bed Frame
Siemens	Mammomat Inspiration
	Siemens RapidPoint Systems
Toshiba	Ultrasound Equipment
Tough-PAC	iPad Covers
Universal Ferrule Remover	Ferrule Remover
Zoll	AED

*Previously Clinell Sporidical Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 16.

www.gamahealthcare.com/compatibility